
DUOWEDGE

www.kasios.com

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
Wedge for high tibial osteotomy (synthetic)

developed in cooperation with Professor Philippe Hernigou. Hôpital Henri Mondor, Créteil - France.


Based on ultra-innovative reconstructive tridimensional ceramic stereolithography technology, DUOWEDGE* is a unique tibial osteotomy wedge featuring two portions with different porosities:

- A solid highly resistant portion, (+80 MPa), which fits into the cortical area, i.e. the maximum load-bearing zone.
- A porous portion to be incorporated in the cancellous area. Osseointegration of this porous section serves to anchor DUOWEDGE in the tibia.

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
- Phosphate tricalcique (TCP) 40%
- Hydroxyapatite (HA) 60%

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
Material

DUOWEDGE is made of phosphocalcic ceramic composed of:

- Tri Calcium Phosphate (TCP) 40%
- Hydroxyapatite (HA) 60%

Porosity

Average porosity of the porous section is approximately 60% with pore size ranging from 400 µm to 700 µm.
The porous network is totally interconnected. It is completely tridimensional and there are no blind pores.

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
Compressive strength*

The compressive strength of the DUOWEDGE solid portion is greater than 80MPa, i.e. approximately 800kg/cm²!
This compressive strength provides reinforcement for the osteosynthesis equipment (plates and screws).

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Test Rm compression
Hernigou

COMPRESSION:
Sensor: 50kN
speed: 0.3mm/min

MAX FORCE: 37500 N
length: 9.930mm
width: 19.91mm
height: 15.01mm

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Packaging

DUOWEDGE units are packaged individually in double sterile blisters.
DUOWEDGE is available in heights of 4mm to 20mm, in increments of 1 mm.

Features

Dual porosity.
Compressive strength greater than 80 MPa.
Osseointegration achieved by the porous portion
Calcium phosphate (identical to mineral phase of bone).
TCP enhances osseointegration process.
Multiple sizes available for more precise valgization.
Rounded edges to enhance ergonomic effect.
Totally interconnected porosity.
Sterile, ready for use.
100% synthetic: no risk of cross-contamination.

Indications

Valgization medial opening wedge osteotomy:

- For medial femorotibial gonarthrosis with genu varum
- For open osteotomy with chronic laxity.
- For open osteotomy to correct genu recurvatum.

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

***Attention!**

DUOWEDGE must always be combined with the use of osteosynthesis equipment.

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Designation	A	B	H	h	Réf./Ref.
Cale/Wedge 4 mm	20	15	4	2	K43504W
Cale/Wedge 5 mm	20	15	5	3	K43505W
Cale/Wedge 6 mm	20	15	6	4	K43506W
Cale/Wedge 7 mm	20	15	7	5	K43507W
Cale/Wedge 8 mm	20	15	8	6	K43508W
Cale/Wedge 9 mm	20	15	9	7	K43509W
Cale/Wedge 10 mm	20	15	10	8	K43510W
Cale/Wedge 11 mm	20	15	11	8,5	K43511W
Cale/Wedge 12 mm	20	15	12	9	K43512W
Cale/Wedge 13 mm	20	15	13	9,5	K43513W
Cale/Wedge 14 mm	20	15	14	10	K43514W
Cale/Wedge 15 mm	20	15	15	11	K43515W
Cale/Wedge 16 mm	20	15	16	11,5	K43516W
Cale/Wedge 17 mm	20	15	17	12,5	K43517W
Cale/Wedge 18 mm	20	15	18	13	K43518W
Cale/Wedge 19 mm	20	15	19	13,5	K43519W
Cale/Wedge 20 mm	20	15	20	14	K43520W

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