kasios®

# JECTOS®

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### Presentation

JectOS is an injectable synthetic bone substitute made of 99% calcium phosphate.

JectOS is a product easy to set and easy to use by surgeons thanks to its long working time (4 to 5 minutes).

JectOS is biocompatible, bioresorbable and osseoconductive. This will ensure proper osseointegration in human bone.

JectOS is delivered in two vials (one vial containing liquid, the other containing powder) and the ancillary accessories necessary for mixing and injecting.

### Indications

JectOS is indicated for the filling of cancellous bone defects.

### Chemical composition

Mixing the powder and the liquid produces a solid compound made of DCPD (55%) and beta TCP (45%). This reaction is slightly exothermic and gives rise to a transient temperature (less than 45°C).

Following the chemical reaction, JectOS is made of two types of calcium phosphate:

- 55% DCPD (Dicalcium phosphate dihydrated)
- 45% TCP (Tricalcium phosphate)

### Physical properties

Porosity: 40%

Pore size: less than 5 µm

Setting time: 9 to 11 minutes (including a working time of 4 to 5 minutes)

Reaction temperature: 43°C – 45°C

Compressive strength: approximately 35 MPa

### Bioresorption

When JectOS is implanted inside a bone defect in close apposition with human bone, it will resorb with time and will be replaced simultaneously by new bone. This affords JectOS a very good osseointegration thanks to its bioresorption and osseoconductivity abilities.

### Features and benefits

- Synthetic: no risk of disease transmission.
- Biocompatible, bioresorbable, bioconductive: good osseointegration.
- Radiopaque: easy to control during injection.
- Injectable: minimally invasive.
- Easy to use.
- Easy to handle: up to 5 minutes working time.
- Putty allows complete filling of the bone defect.
- Sold sterile.
- Slight temperature rise: between 43°C and 45°C.

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
**1. Pour entire quantity of liquid into bowl.**  
Caution! Be sure to empty the complete vial. Do not leave any liquid inside.

**2. Add entire quantity of powder contained in the other vial.**  
Caution! Be sure to empty the complete vial. Do not to leave any powder inside.

### START THE CLOCK!

**3. Combine powder and liquid together and mix vigorously for 30 seconds until mixture becomes a smooth homogeneous liquid.**

**4. Aspirate the liquid through the cannula connected to the syringe.**

**5. From time to time, test a drop of the product to check for hardening.**

**6. When JectOS reaches a toothpaste like consistency, start injecting into the bone defect.**

**CAUTION!**  
Do not implant JectOS before it looks like toothpaste!

As the chemical reaction is highly temperature dependent, setting time can vary. That is why it is very important to check the hardening of the product from time to time between minute 1 and minute 5.

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

◇ **73-year-old woman**  
**Fracture of the distal radius.**

Cominución circunferencial.  
Radiografías preoperatorias y foto peroperatoria.  
Circumferential comminution.  
Pre-operative X-rays and peri-operative palmar view.

Osteosíntesis con placa.  
Radiografía postoperatoria inmediata.  
Relleno de la cominución circunferencial (1).  
Atornillamiento en el cemento (2).  
Osteosynthesis using a palmar plate.  
Immediate post-operative X-rays.  
Filling the circumferential comminution (1).  
Screwing into the cement (2).

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
◇ **26-year-old man**  
**Extra-articular malunion at one year.**

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◇ 64-year-old man  
Filling the distal radius where can-  
cellous bone had been removed  
(thumb reconstruction).

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◇ 58-year-old man  
Extra articular malunion at 6 months.  
Fracture initially by a simple extra fixator.

Osteosíntesis inicial realizada únicamente con fijador externo.  
Fracture initially stabilized by a simple external fixator.

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Each box contains:

- 1 vial of powder
- 1 vial of liquid
- 1 mixing bowl
- 1 spatula
- 1 syringe
- 1 cannula

DISTRIBUIDO POR  
DISTRIBUTED BY

made by  
**kasios®**  
"Specializing in bone substitutes"

18 CHEMIN DE LA VIOLETTE - 31240 L'UNION - FRANCE  
Tél. +33 (0)5 34 27 33 23 . Fax +33 (0)5 34 27 33 27  
contact@kasios.com

**www.kasios.com**